

## 2019 Novel Coronavirus

# Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens from Persons Under Investigation (PUIs) for 2019 Novel Coronavirus (2019-nCoV)

February 2, 2020

Health care providers should contact their local/state health department immediately to notify them of patients with fever and lower respiratory illness who traveled to Wuhan, China within 14 days of symptom onset. Local and state public health staff will determine if the patient meets the [criteria for a person under investigation \(PUI\)](#) for 2019 Novel Coronavirus (2019-nCoV). Clinical specimens should be collected from PUIs for routine testing of respiratory pathogens at either clinical or public health labs. Note that clinical laboratories should NOT attempt viral isolation from specimens collected from 2019-nCoV PUIs.

At this time, diagnostic testing for 2019-nCoV can be conducted only at CDC.

State and local health departments who have identified a PUI should immediately notify CDC's Emergency Operations Center (EOC) at 770-488-7100 to report the PUI and determine whether testing for 2019-nCoV at CDC is indicated. The EOC will assist local/state health departments to collect, store, and ship specimens appropriately to CDC, including during afterhours or on weekends/holidays.

Testing for other respiratory pathogens by the provider should be done as part of the initial evaluation and should not delay specimen shipping to CDC.

If a PUI tests positive for another respiratory pathogen, after clinical evaluation and consultation with public health authorities, they may no longer be considered a PUI.

## Specimen Type and Priority

For initial diagnostic testing for 2019-nCoV, CDC recommends collecting and testing upper respiratory (nasopharyngeal AND oropharyngeal swabs), and lower respiratory (sputum, if possible)) for those patients with productive coughs. Induction of sputum is not indicated. Specimens should be collected as soon as possible once a PUI is identified, regardless of the time of symptom onset. Maintain [proper infection control](#) when collecting specimens.

## General Guidelines

Store specimens at 2-8°C and ship overnight to CDC on ice pack. Label each specimen container with the patient's ID number (e.g., medical record number), unique specimen ID (e.g., laboratory requisition number), specimen type (e.g., serum) and the date the sample was collected. Complete a [CDC Form 50.34](#) for each specimen submitted. In the upper left box of the form, 1) for *test requested* select "Respiratory virus molecular detection (non-influenza) CDC-10401" and 2) for *At CDC, bring to the attention of* enter "Stephen Lindstrom: 2019-nCoV PUI".

## I. Respiratory Specimens

### A. Lower respiratory tract

#### Bronchoalveolar lavage, tracheal aspirate

Collect 2-3 mL into a sterile, leak-proof, screw-cap sputum collection cup or sterile dry container. Refrigerate specimen at 2-8°C and ship overnight to CDC on ice pack.

#### Sputum

Have the patient rinse the mouth with water and then expectorate deep cough sputum directly into a sterile, leak-proof, screw-cap sputum collection cup or sterile dry container. Refrigerate specimen at 2-8°C and ship overnight to CDC on ice pack.

## B. Upper respiratory tract

### Nasopharyngeal swab AND oropharyngeal swab (NP/OP swab)

Use only synthetic fiber swabs with plastic shafts. Do not use calcium alginate swabs or swabs with wooden shafts, as they may contain substances that inactivate some viruses and inhibit PCR testing. Place swabs immediately into sterile tubes containing 2-3 ml of viral transport media. NP and OP specimens should be kept in separate vials. Refrigerate specimen at 2-8°C and ship overnight to CDC on ice pack.

*Nasopharyngeal swab:* Insert a swab into the nostril parallel to the palate. Leave the swab in place for a few seconds to absorb secretions. Swab both nasopharyngeal areas with the same swab.

*Oropharyngeal swab (e.g., throat swab):* Swab the posterior pharynx, avoiding the tongue.

### Nasopharyngeal wash/aspirate or nasal aspirate

Collect 2-3 mL into a sterile, leak-proof, screw-cap sputum collection cup or sterile dry container. Refrigerate specimen at 2-8°C and ship overnight to CDC on ice pack.

## II. Serum

**Minimum volume required:**

*Children and adults:* Collect 1 tube (5-10 mL) of whole blood in a serum separator tube.

*Infant:* A minimum of 1 mL of whole blood is needed for testing pediatric patients. If possible, collect 1 mL in a serum separator tube.

Serum separator tubes should be stored upright for at least 30 minutes, and then centrifuged at 1000–1300 relative centrifugal force (RCF) for 10 minutes before removing the serum and placing it in a separate sterile tube for shipping (such as a cryovial). Refrigerate the serum specimen at 2-8°C and ship overnight to CDC on ice-pack.

## III. Shipping

Specimens PUI's must be packaged, shipped, and transported according to the current edition of the [International Air Transport Association \(IATA\) Dangerous Goods Regulations](#) [↗](#). Store specimens at 2-8°C and ship overnight to CDC on ice pack. If a specimen is frozen at -70°C ship overnight to CDC on dry ice. Additional useful and detailed information on packing, shipping, and transporting specimens can be found at [Interim Laboratory Biosafety Guidelines for Handling and Processing Specimens Associated with 2019 Novel Coronavirus \(2019-nCoV\)](#).

*For additional information, consultation, or the CDC shipping address, contact the CDC Emergency Operations Center (EOC) at 770-488-7100.*